# Geratherm<sup>®</sup> chlamydia check

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### Produktbeschreibung

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
### Product description

#### Delivery scope / Package contents

- Aluminum bag with test cassette and desiccant
  - Buffer bottle with buffer solution 0,6 ml
  - Swab (manufacturer: Puritan Medical Products Company LLC), individually sterile packed
  - Instructions for use
- You will also need a watch.

### Product overview / description of the components

### Mode of operation / measuring principle

The Geratherm<sup>®</sup> chlamydia check for self-diagnosis is an immunological, qualitative lateral flow rapid test for the detection of Chlamydia trachomatis (antigen) in cervical swab samples. When the sample (vaginal secretion mixed with buffer solution) is applied to the sample field „S“, it migrates along the membrane by capillary action. In a positive sample, an anti-chlamydia trachomatis antibody-gold complex binds the chlamydia trachomatis antigen. This complex binds to an immobilized anti-chlamydia-trachomatis antibody on the membrane at the level of the „T“ mark of the test cassette. This results in the aggregation of the antibody-gold complex, producing a purple line in the result field at „T“. A colored line at the level of the control field „C“ serves as a procedural control, indicating that the correct amount of sample was applied and that membrane wetting was sufficient.

### Clinical benefit

This test serves to detect Chlamydia trachomatis, one of the three known types of Chlamydia. Chlamydia infection is the most frequently occurring sexually transmitted disease. In Germany alone, the number of new infections is estimated at between 300,000 and 500,000 per year. Younger adults (between age 15 and 25), those who frequently change sexual partners, the sexual partners of those infected (even those without symptoms who are unaware that they carry the infection), as well as newborn babies of infected mothers, are the most frequently affected. Infections without sexual contact (e.g. from swimming pools) are rare, but not impossible. Chlamydia trachomatis infections frequently have no symptoms. However, they may also lead to complaints such as eye infections, arthritis, and urinary-tract infections. In addition, untreated infection frequently leads to infertility, as well as increasing the risk of ectopic pregnancy and preterm births. In developing countries, Chlamydia trachomatis infection is the most common cause of blindness among newborns. It is particularly important to treat Chlamydia infections in pregnancy. Around 50 % of infected mothers pass the infection on to their newborns. Common complications among newborns include conjunctivitis and pneumonia.

### Active ingredients

- Monoclonal mouse anti-chlamydia tr. Antibody colloidal gold conjugate (3,2 µg/ Test) + Polyclonal rabbit IgG antibody colloidal gold conjugate (1,94 µg/Test)
- Monoclonal mouse anti-chlamydia tr. Antibody (2,6 µg/Test)
- Polyclonal goat anti-rabbit IgG Antibody (5 µg/Test)

### Performance data

This test achieved a diagnostic sensitivity of 85,7 % and a diagnostic specificity of 98,3 % in the clinical performance evaluation with 596 participants compared to the PCR method. There are no known cross-reactions. This test detects the 15 serotypes of the Chlamydia trachomatis bacterium. The detection limit for this test is 2,45 µg/ml elemental body. In addition, a lay study was conducted with 56 subjects. Of these 56 subjects, 54 (96,4 %) correctly interpreted the result.

### Cross-reactions:

No cross-reactions are expected with the microorganisms tested at the following concentrations:

Concentration: 10 <sup>7</sup> CFU/ml	Aeromonas spp; Campylobacter spp; Citrobacter spp; Enterobacter spp; Haemophilus coli; Herpes simplex spp; Lactobacillus spp; Mycoplasma spp; Neisseria meningitidis; Proteus spp; Salmonella spp; Shigella spp; Trichomonas spp; Veillonella spp; Bacteroides spp; Candida spp; Clostridium spp; Escherichia coli; Gardnerella spp; Klebsiella spp; Listeria spp; Neisseria gonorrhoeae; Peptococcus spp; Pseudomonas spp; Serratia spp; Streptococcus spp; Ureaplasma urealyticum; Yersinia spp
Concentration: 10 <sup>9</sup> CFU/ml	Staphylococcus spp (coag. neg.); Staphylococcus spp (coag. pos.)

### Interference Substances:

No interference with the test is expected with the tested substances.

### Potentially interfering substance

Capillary blood (~3,3 %)  
Seminal fluid  
Mucus  
Haemorrhoidal creme (Hädensa<sup>®</sup>)  
Leukozytes  
Vaginal medical products: Canesten<sup>®</sup> (Clotrimazol), Lomexin<sup>®</sup> (Fenticonazol)  
Vaginal suppositories: Multi-Gyn Actigel<sup>®</sup> (Galactoarabinan, Polyglucuronid Acid)

### Intended purpose

Self-test to detect the bacterium chlamydia trachomatis

### Indications

- Chlamydia trachomatis is the most common sexually transmitted pathogen. The Geratherm<sup>®</sup> chlamydia check should therefore be used regularly by women,
- if a sexually transmitted infection is suspected
- in advance to pregnancy
- who have sex with new or frequently changing partners
- who are sexually active.

Furthermore, the Geratherm<sup>®</sup> chlamydia check is indicated for women who experience the following symptoms:

- pain when urinating
- unusual vaginal discharge
- pain in the tummy or pelvis
- pain during sex
- bleeding after sex
- bleeding between periods

### Important note

This test is not suitable for use during pregnancy or during menstruation or for three days after the end of your period or if you have a urinary tract infection.

### Target patient group(s) and intended users

This test is intended for sampling from the cervix, i.e., exclusively for self-testing in women.

### Notes, precautions and warnings

- Please read the instructions for use in detail before using the test.
- These instructions must always be followed.
- Do not use after the expiration date indicated on the aluminum bag.
- Only open the aluminum bag before performing the test.
- Do not use if the aluminum bag is open or damaged.
- Do not use if components are damaged.
- For in vitro use only.
- Use test outside the body only.
- Keep away from children.
- Use only the enclosed swab and make sure that the original packaging is undamaged, or sterility is not compromised before using it!
- The buffer solution in the bottle is corrosive. Avoid skin contact!
- There is deliberately little buffer solution in the buffer bottle.
- Do not damage or disassemble the test cassette!
- Greater amounts of haemoglobin in the sample can lead to false positive results.
- To be safe, wash your hands after performing the test.
- Do not make any medically important decisions without your physician!
- This test is subject to a degree of uncertainty. A physician should always be consulted.
- This test is for sampling from the cervix (neck of the womb), i.e. for use by women only.
- The test materials of animal origin used (e.g. antibodies) are potentially infectious materials, but there is no danger from them if you use all test components according to the instructions for use.
- Only use components that are included in the packaging and make sure that the original packaging is undamaged before use!
- Do not use this test if you have any injuries in the sampling area.
- This test is intended for single use only.
- Do not use this test after the expiration date.
- When using the buffer solution, apply only 2 drops of the sample solution to the sample field.
- Do not exceed the waiting time to evaluate the test.
- Do not use the test during pregnancy.
- In case of infection, it is recommended that the sexual partner also undergoes testing and treatment.

### Use

#### Preparation / Taking a sample using the swab

Before taking a sample, avoid urinating for at least one hour, since otherwise an incorrect result may be yielded. Perform the test in a location where you will be undisturbed. Please put the aluminum bag aside UNOPENED for the time being, this will only be opened and needed to carry out the test. Take the swab out of the pack. Only use the provided swab and make sure that the original packaging is undamaged before using it! If the packaging is damaged, please do not use the swab and instead contact the manufacturer.


- Select a comfortable position, such as one you would adopt for inserting a tampon. Carefully insert the swab into your vagina as far as the neck of the womb (top end of your vagina). Rotate the swab evenly for about 30 seconds and then slowly withdraw the swab! Wipe off any excess vaginal mucus before taking the sample!

- There is purposely only a small quantity of solution in the sample bottle. Dip the swab containing the sample right into the bottle so that it touches the bottom of the sample bottle. Rotate the swab in the solution and squeeze it out several times against the side wall of the bottle in order to ensure that the sample is mixed thoroughly with the buffer solution.

- Wait 5 minutes. If you wait too long, or not long enough, this will distort the test result. During the waiting time, the swab must remain in the bottle.

### Test execution

- Now prepare the test cassette. To do this, open the aluminum bag, remove the test cassette and place it on a clean flat surface. In only 15 minutes you will find out whether an infection with Chlamydia trachomatis is present.

- Now wait 15 minutes.** If you wait for longer than 20 minutes, it will no longer be possible to interpret the result reliably.

- To be on the safe side, wash your hands after performing the test.

### Reading and interpreting the result

A sample is considered **POSITIVE** if, within the reaction time of 15 minutes, two purple lines appear in the reaction zone of the cassette (at „C“ and „T“), even if the coloured line at „T“ is very faint. A positive result can often be identified before the end of the 15-minute reaction time.

A result is considered **NEGATIVE** if there is only one purple line at „C“. A negative result can only be established with certainty at the very end of the 15-minute reaction time.

A result is considered **INVALID** if no purple line appears or there is only a line at „T“ in the test field, or the entire test field turns purple.

Possible reasons for an invalid result may be a damaged foil or swab pack, incorrect storage, or an error in performing the test. Please retain all the test materials and contact the manufacturer.

**If you wait more than 20 minutes, it will no longer be possible to evaluate the test reliably. Please therefore dispose of the test cassette once you have finished interpreting the result at the correct time.**

### Question and answers (Disclosure of residual risk)

#### What should I do if my test result is positive?

Please see your doctor in order to discuss the result. Please take these instructions for use into your appointment, so that your doctor is able to gain a clearer picture of the test you have performed.

#### What should I do if my test result is negative?

A negative test result means that you are highly unlikely to have a Chlamydia trachomatis infection, but it is still recommended that you are regularly checked by your doctor.

### Literature

- Cristillo et al. Point-of-Care Sexually Transmitted Infection Diagnostics: Proceedings of the STAR Sexually Transmitted Infection—Clinical Trial Group Programmatic Meeting. Sexually Transmitted Diseases. 2017 - Volume 44 - Issue 4 - p 211-216.
- Langille et al. A pilot project for chlamydia screening in adolescent females using self-testing: characteristics of participants and non-participants. Can J Public Health. 2008 Mar-Apr;99(2):117-20.
- Greer et al. Rapid Diagnostic Methods in Sexually Transmitted Infections. Infect Dis Clin N Am, 22 (2008) 601–617.
- Rodrigues et al. Chlamydia trachomatis as a Current Health Problem: Challenges and Opportunities. Diagnostics (Basel). 2022 Jul 25;12(8):1795.

### Storage and Transport

Store in a dry place between 4 °C and 28 °C.

### Disposal of packaging and consumables

All test components (including the desiccant) can be disposed of in household waste after performance. Make sure that individual packaging components are disposed of in accordance with local regulations.

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




- Cristillo et al. Point-of-Care Sexually Transmitted Infection Diagnostics: Proceedings of the STAR Sexually Transmitted Infection—Clinical Trial Group Programmatic Meeting. Sexually Transmitted Diseases: 2017 - Volume 4 - Issue 4 - p 211-218.
- Langille et al. A pilot project for chlamydia screening in adolescent females using self-testing: characteristics of participants and non-participants. Can J Public Health. 2008 Mar-Apr;99(2):117-20.
- Greer et al. Rapid Diagnostic Methods in Sexually Transmitted Infections. Infect Dis Clin N Am. 22 (2008) 601–617.
- Rodrigues et al. Chlamydia trachomatis as a Current Health Problem: Challenges and Opportunities. Diagnostics (Basel). 2022 Jul 25;12(8):1795.

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## CE0197

	Hisopo aplicador Puritan Medical Products Company LLC 31 School Street Guilford EE. UU.		
			

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