Name / Name Nom / Nome		Sandra Mustermann		Monat / Month Mois / Mese		aktueller Monat / current month / Mois en cours Mese corrente	
Zykluslänge / Cycle length Durée du cycle / Durata del ciclo		28 Tage / Days / Jours / Giorni		Normale Länge der Periode / Usual period length / Durée habituelle des règles / Normale durata delle mestruazioni		6 Days	
Messort / Measurement site Site de mesure / Sede di misurazione		oral / oral buccal / orale		rektal / rectal rectal / rettale		vaginal / vaginal vaginal / vaginale	
Kalendertag / Calendar day / Jour du calendrier / Giorno di calendario		18 19 20 21 22 23 24 25 26 27 28 29 30 31 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40					
Zyklusstag / Cycle day / Jour du cycle / Giorno del ciclo		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40					
		37,50					
		37,40					
		37,30					
		37,20					
		37,10					
Basal- bzw. Aufwachtemperatur in °C		37,00					
		36,90					
Basal or waking temperature in °C		36,80					
		36,70					
Température basale ou au réveil en °C		36,60					
		36,50					
Temperatura basale o al risveglio in °C		36,40					
		36,30					
		36,20					
		36,10					
		36,00					
		35,90					
Menstruation / Menstruation Menstruation / Mestruazione		XXXXXX					
Geschlechtsverkehr / Intercourse Rapport sexuel / Rapporti sessuali				X XXXX			
Zervixschleim / Cervical mucus Glaire cervicale / Muco cervicale				XXX			
Sonstiges / Other Autre / Altro							

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
FOR FAMILY PLANNING & CYCLE TRACKING

The Geratherm® *basal* supports you in natural family planning and cycle tracking. By measuring your basal body temperature every day, your individual temperature pattern gradually becomes visible. This helps you reliably identify ovulation and your fertile days. After about three cycles, your personal temperature pattern is usually clearly recognizable.

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HOW TO MEASURE CORRECTLY

1. Measure immediately after waking up, before getting out of bed, and preferably at the same time every morning.
2. Before measuring, shake the thermometer down with the reset aid until the liquid column is below 36 °C.
3. Choose one measurement site (oral, vaginal or rectal) and use the same site throughout the entire cycle.
4. Measure for the same length of time every morning, at least four minutes. Only a consistent measuring time ensures comparable results.
5. Read and record the temperature. The reading remains visible until the thermometer is reset, so you can record it at your convenience.

Tip: Sleep for at least five hours and avoid eating, drinking or smoking before measuring. Clean the thermometer before and after each use.

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