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Medical digital thermometer Geratherm® Color GT-131

Geratherm Medical GERMANY GmbH Germany CE 0197 

PARTS DESCRIPTION

INSTRUCTION FOR USE

Please read carefully before using


SPECIFICATIONS

Range: 32,0 °C to 43,9 °C, Below 32,0 °C display L °C, Above 43,9 °C display H °C
 Intended use: Measurement of body temperature
 Resolution: 0,1 °C
 Accuracy: ± 0,1 °C (34,0 °C to 42,0 °C) ± 0,2 °C (other range)
 Display: Liquid crystal display 3 1/2 digits
 Battery (included): Micro Alkaline 192 or LR41 1,55V
 Power consumption: 0,15 milliwatts in measurement mode
 Battery life: More than 200 hours of continuous operation
 IP classification: IP 27 (Protected against solid foreign objects with a diameter of 12,5 mm and protection against temporary immersion)
 Dimensions: 129,80 mm x 18,00 mm x 10,00 mm
 Weight: Approx. 14 grams including battery
 Beeper: Approx. 10 secs. sound signal when no relevant increase in temperature is detected
 Memory: Built-in memory enabling display of the previously measured value
 Working Conditions: Temperature: 10 °C to 40 °C, Humidity: 15 % to 95 % R.H.
 Storage Conditions: ⚠ It is only possible to measure temperatures that are above the ambient temperature. Temperature: -25 °C to +55 °C, Humidity: 15 % to 95 % R.H.
 Guarantee of Quality: The Geratherm® color complies requirements DIN EN ISO 80601-2-56: Medical electrical equipment. Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement (ISO 80601-2-56). The Geratherm Medical GERMANY GmbH is certified according to DIN EN ISO 13485 and complies with Council Directive 93/42/EEC. This entitled to affix the CE 0197 marking (Notified body: TÜV Rheinland LGA Products GmbH). The present medical product is a device with extremely high degree of noise immunity. Therefore, the enclosure with the EMC tables was not printed out.

PRECAUTIONS

- Use of this instrument is not intended as a substitute for consultation with your physician. Prolonged high fever requires medical attention especially for your children. Please contact your physician.
- Child's temperature shall be taken under adult's surveillance.
- Keep out of the reach of unsupervised children. Choking resulting from a child swallowing a small part that has become detached from the thermometer.
- Do not boil, bite, bend, drop or disassemble this thermometer.
- Do not attempt to disassemble the unit except to replace the battery.
- Avoid taking temperature in 30 minutes after physical exercises, bathing or dining.
- Do not store, use and transport the unit under direct sunlight, at a high temperature, or in high humidity or dust. Performance may be degraded.
- Stop using the thermometer if it operates erratically or if the display malfunctions.
- If the thermometer has been stored at below-freezing temperatures, allow the thermometer to warm naturally to room temperature before using.

BATTERY REPLACEMENT

- When the symbol  appears, the battery needs to be replaced.
- Pull the battery cover and the plug off in the direction as shown.
- Use a toothpick or insulated probe to remove the battery. Avoid using any sharp metal object in this operation.
- Place a new 1,55 VDC button size battery type SR41, UCC392, LR41, 192 or equivalent in the chamber with positive side up and the negative side facing down.
- Replace and fasten the plug and battery cover securely.

⚠ Discard old battery carefully, out of reach of young children. Do not dispose of the battery in a fire. It may explode.

HOW TO USE

- Press the on/off button. The display will read .
- Release the power switch and the display will show L °C flashing.
- Clean or disinfect the probe before every measurement. Position the probe (oral, rectal, axillary etc.).
- As soon as the thermometer does not register any further relevant increase in temperature, °C stops flashing and the maximum value measured up to that time is shown, accompanied by a series of audible signals.
- The thermometer continues measuring the temperature for a further minute after the audible signals begin, and in this time there may possibly be a further slight increase in the temperature shown. After that, the maximum temperature measured is recorded and stored; it remains visible on the LCD display for 3 minutes until the thermometer is automatically switched off. Or short press the ON/OFF button to switch off directly.
- The unit will automatically turn off in 3 minutes (approx.). However, to prolong battery life, it is best to turn the thermometer off by pressing the ON/OFF button once the temperature has been noted.
- To check memory: Press on/off button and hold button for 2 seconds when the unit is under off state. The last temperature will be recalled on the LCD display.
- NOTE:** You may put a sanitary plastic probe shield on the sensor tip to help prevent cross-infection and recontamination. Use of the probe cover may result in a 0,1 °C difference from actual temperature. Dispose of any used probe shield and store the unit in a safe place.

MEASUREMENT METHODS

It is only possible to measure temperatures that are above the ambient temperature.

ORAL USE (about 45 – 60 seconds)

Place the probe well under the patient's tongue. Instruct the patient to keep their mouth closed. A normal temperature by this method is usually considered between 36.0 °C and 37.5 °C.

RECTAL USE (about 30 – 60 seconds)

If rectal temperature is recommended by your doctor, you may use a STOP cover. Insert thermometer into the probe cover. Lubricate the probe cover with a water soluble jelly for easier insertion. Do not use petroleum jelly. Insert tip of probe no more than 1,5 cm into rectum. STOP if you meet any resistance. In general, the rectal temperature is around 0.5 °C higher than oral temperature.

AXILLAR USE

Wipe armpit with a dry towel. Place probe in the patient's armpit and keep the patient's arm pressed firmly against their side. In general, the axillary temperature is around 0.5 °C lower than oral temperature. Note that in order to achieve better axillary temperature measurement, a minimum measuring time of 2 to 3 minutes is recommended regardless of the beep sound.

If you compare temperatures between two thermometers, it is possible to have two different readings even if both thermometers are within their accuracy tolerances. Differences in how you take each reading and the thermometer's tolerance range can account for nearly 0.1 °C variability.

REMEMBER: Find your normal temperature and watch for changes.

CLEANING AND DISINFECTION

For household use, always clean the probe head before and after each use. Wash the tip with a solution of mild detergent and cool water. Disinfect the thermometer by wiping the sensor and lower stem with a cloth dipped in a household antiseptic solution such as rubbing alcohol (recommendation: 70% isopropyl alcohol). The thermometer is water resistant. Do not boil or clean in autoclave.

DISPOSAL

Observe the applicable regulations when disposing of the thermometer and batteries. This thermometer must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the thermometer. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

LIMITED WARRANTY

This thermometer is guaranteed for two years from the date of purchase against manufacturer's defect under normal use. If your unit does not function properly due to defective parts or assembly we will repair it free of charge. All parts are covered by this warranty excluding the battery. The warranty does not cover damage to your unit due to improper handling. We recommend that the accuracy of the unit be checked every 2 years by an authorised laboratory.

SYMBOL INDEX

	Follow instructions for use		Type BF applied part		Reference number		Model number
	Keep dry		Manufacturer		Batch code (mm/yyyy; month/year)		Medical device
	Store between -25 °C and +56 °C		95 % Store between 15 % and 95 % R.H.		The device must not be disposed of with household waste.		Voltage of battery

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