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**BIBLIOGRAPHIE**

- Henry JB, et al. Clinical Diagnosis and Management by Laboratory Methods, 20th Ed. Philadelphia, Saunders. 371-372, 375, 379, 382, 385, 2001.
- Ma Junlong, Cong Yulong. The effect of bacteriuria on the determination of urine red blood cells by urine analyzer. Chinese Journal of Medical Examination, 1999, 22(4): 205.
- Shuai Lihua, Jiujiang Medical Journal 2002, 17 (2): 122.

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18

**MATERIALS PROVIDED**

Test strip, Plastic cup, Color chart, Package insert

**MATERIALS REQUIRED BUT NOT PROVIDED**

Timer or a watch with a second hand, Specimen container

**PROCEDURE**

**ATTENTION: It's recommended to take a sample of urine for the test in the early morning since it's the most concentrated.** The urine used for the test should not come into contact with water from the toilet or any disinfectant or cleaning substances.

**For women only: The test should not be performed during or for three days after your menstrual period.** The urine sample should not be contaminated with vaginal fluids since this may produce a misleading result.

**Do not make any important medical decision without first referring to your doctor.**

**COLLECT URINE:** Collect part of the urine in the supplied plastic cup or using a clean cup without any residual detergents. Make sure to fill up the cup with urine.

**PERFORMING THE TEST:**

1) Open the foil pouch and take out the test strip. **Do not touch the test fields.** Once opened the pouch, it is recommended to perform the test immediately.

2) Dip the test strip in the urine sample. **ATTENTION:** Press the strip and make sure that all four test fields are immersed for about **1-2 seconds**.

3) Then remove the test strip and wipe off any surplus urine against the rim of the container or with an absorbent material (e.g. a paper towel) to avoid mixing chemicals from adjacent reagent areas.

4) **Wait for 2 min.** (do not read results after 3 minutes) **Read the result separately for each parameter, compare color with color chart provided.**

**READING THE RESULTS**

**Read the result separately for each parameter, compare color with color chart provided.** Color changes on the edges of the test pads or color changes after more than 3 minutes have to be ignored.

**NEGATIVE**  
The test field for **LEUKOZYTEN** stayed whitish. The test field for **BLOOD** stayed mustard yellow.

**POSITIVE RESULT FOR LEUKOZYTEN**  
If the color of the test field has changed to purple, then leukocytes have been found in your urine.

**POSITIVE RESULT FOR BLOOD**  
If the color of the test field has changed to green, then some green spots appear on a white background, then blood has been found in your urine.

**POSITIVE RESULT FOR NITRITE**  
If the color of the test field has changed to pink, then nitrites have been found in your urine.

**POSITIVE RESULT FOR PROTEIN**  
If the color of the test field has changed to green, then proteins have been found in your urine.

**TECHNICAL NOTES ON PARAMETERS**

The Test detects **LEUKOZYTEN, BLOOD, NITRITE** and/or **PROTEINE** in urine.

**LEUKOZYTEN:** The presence of leukocytes in urine is an important symptom of an inflammation of the kidneys and the urinary tract, protein react with the pad and changing its color to purple.

When taking cephalaxin and cephalothin, or high concentration of oxalic acid may also cause test results to be artificially low. Tetracycline may cause decreased reactivity, and high concentrations of oxalic acid may also cause test results to be artificially low.

**BLOOD:** A uniform green color conversion indicates the presence of hemoglobin or hemolyzed erythrocytes; scattered or compacted green spots indicate intact erythrocytes. Generally urine occult blood attributed to the following three reasons, one is stone, inflammation, cancer. On the aspect of inflammation, such as glomerulonephritis, pyelonephritis, cystitis, but may have menstruation. urine occult blood will have occurred. Whether stones, kidney, ureter or bladder stones, may cause other situations such as occult blood. The tumor can also cause occult blood, such as benign or malignant tumor of kidney, ureter and bladder.

Menstrual period, constipation may cause a positive result.

**NITRITE:** Gram-negative bacteria in urine convert nitrate from food into nitrite. Nitrite reacts with a chemical in the test field and leaves a pink shade. The test result may be distorted if urine does not stay for long in the bladder, due to hunger, a vegetable-free diet or antibiotic treatment. Comparing the test on a white background may aid in the detection of low nitrite levels, which might otherwise be missed.

**PROTEIN:** An indicator on the test field reacts with protein in the urine, changing its color to green. They may be found where there is inflammation of the bladder or prostate or bleeding in the urinary tract. Infusions containing polyvinylpyrrolidone may yield a false positive result. Chemical components in the test fields must be viewed as potentially dangerous substances, although they present no hazard provided that all test components are used in accordance with these instructions.

**CONTROL PROCEDURE**

To serve as good quality control, the instruction must be followed closely when performing the test. Failure to follow directions in insert may yield inaccurate test results.

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Geratherm Medical GERMANY GmbH Fahrenheitsstraße 1 99331 Geratal/Germany